



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse COX411 Monoclonal Antibody

单克隆抗体

产品货号: MA01931

第三版

免疫原:Fusion protein of COXIV

亚型:IgG1

克隆:550

描述:COX411,also named as COX4 and COXIV-1,belongs to the cytochrome c oxidase IV family. It is one of the nuclear-coded polypeptide chains of cytochrome c oxidase,the terminal oxidase in mitochondrial electron transport. COX411 is a marker for mitochondria. It has two isoforms (isoform 1 and 2). Isoform 1(COX411) is ubiquitously expressed and isoform 2 is highly expressed in lung tissues. COX411 is commonly used as a loading control. This antibody is specific to COX411 and do not cross reacts with COX412.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF