

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CSK Monoclonal Antibody

单克隆抗体

产品货号: MA01935

第三版

免疫原:Fusion protein of CSK

亚型:IgG2b

克隆:240

描述:CSK(Tyrosine-protein kinase CSK) is also named as C-Src kinase,protein-tyrosine kinase CYL and belongs to the protein kinase superfamily. Csk is a non-receptor tyrosine kinase with a molecular mass of 50 kDa. It contains the SH3 and SH2 domains in its N-terminus and a kinase domain in its C-terminus. It negatively regulates the kinase activity of SFKs and plays an important role in physiological functions via signaling pathways for cell proliferation,differentiation,adhesion,and migration. Csk may play a role in tumor suppression by inhibiting the oncogenic activity of SFKs. Indeed,overexpression of Csk can inhibit tumor growth of human colon cancer cells. This protein can be autophosphorylated.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF