

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse c-SRC Monoclonal Antibody

单克隆抗体

产品货号: MA01936

第三版

免疫原:Fusion protein of c-SRC

亚型:IgG2b

克隆:38

描述:SRC,also named as SRC1 and p60-Src,belongs to the protein kinase superfamily,Tyr protein kinase family and SRC subfamily. It is a non-receptor protein tyrosine kinase that plays pivotal roles in numerous cellular processes such as proliferation,migration,and transformation. In concert with PTK2B, SRC plays an important role in osteoclastic bone resorption. It promotes energy production in osteoclasts by activating mitochondrial cytochrome C oxidase.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF