

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CXCR4 Monoclonal Antibody

单克隆抗体

产品货号: MA01940

第三版

免疫原:Fusion protein of CXCR4

亚型:IgM

克隆:226

描述:The chemokine receptor CXCR4,a 7-transmembrane protein on cell surface,acts with the CD4 protein to support HIV entry into cells. CXCR4 is highly expressed in brain and heart,white blood cells,vascular endothelial cells,and umbilical cord endothelial cells. Several CXCR4 isoforms have been identified with molecular weights of 47,50,62 and 98kDa. The 62-64kDa CXCR4 is predominantly ubiquitinated and is preferentially associated with CD4. This antibody recognizes the 62kd CXCR4.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF