



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CYB5R3 Monoclonal Antibody

单克隆抗体

产品货号: MA01942

第三版

免疫原:Fusion protein of CYB5R3

亚型:IgG2b

克隆:

描述:CYB5R3,also named as DIA1 and B5R,belongs to the flavoprotein pyridine nucleotide cytochrome reductase family. It plays an essential role in the mitochondrial metabolism by its contribution to cellular redox homeostasis. Overstimulation of CYB5R3 associated with lipid rafts can account for the overshoot of plasma membrane-focalized superoxide anion production that triggers the entry of CGN in the irreversible phase of apoptosis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF