



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Cytoglobin Monoclonal Antibody

单克隆抗体

产品货号: MA01945

第三版

免疫原:Fusion protein of Cytoglobin

亚型:IgG1

克隆:77

描述:Cytoglobin is a ubiquitously expressed hexacoordinate hemoglobin that may facilitate diffusion of oxygen through tissues. The protein,with calculated MW of 21 kd,runs as a 29kd one in canine retina,kidney,liver,lung,and heart. It may be caused by posttranslational modification of the protein. Various immuno-staining patterns were reported including nuclear,cytoplasm and extracellular location.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF