

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CK-19 Monoclonal Antibody

单克隆抗体

产品货号: MA01947

第三版

免疫原:Fusion protein of Cytokeratin 19

亚型:IgG2b

克隆:335

描述:The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells. KRT19,one of type I keratins,is specifically expressed in the periderm,the transiently superficial layer that envelopes the developing epidermis. Due to its high sensitivity,KRT19 is the most used marker for detection of tumor cells disseminated in lymph nodes,peripheral blood,and bone marrow of breast cancer patients. This antibody,generated against full length KRT19 protein,also recognizes KRT17,another type I keratin homologous to KRT19. KRT17 is used as a marker for trauma.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF