



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse DSC2 Monoclonal Antibody

单克隆抗体

产品货号: MA01952

第三版

免疫原:Fusion protein of Desmocollin 2

亚型:IgG2a

克隆:118

描述:Desmocollin 2 (DSC2) is a member of the desmocollin subfamily of the cadherin superfamily. The desmocollins (DSC1-3), along with the desmogleins (DSG1-4), are found primarily in epithelia and cardiac muscle where they constitute the adhesive proteins of the intercellular desmosome junctions and are required for cell adhesion and desmosome formation. Mutations in DSC2 gene are associated with arrhythmogenic right ventricular dysplasia-11.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF