



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ECM1 Monoclonal Antibody

单克隆抗体

产品货号: MA01956

第三版

免疫原:Fusion protein of ECM1

亚型:IgG2a

克隆:480

描述:Extracellular matrix protein 1 (ECM1) is a glycoprotein involved in a number of biological processes such as bone formation, skin differentiation, cell proliferation, and promotes angiostasis. Pathologically, ECM1 contributes to the formation and metastasis of several types of cancer including breast, thyroid and hepatocellular cancers. (21128013) It inhibits MMP9 proteolytic activity.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF