

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat EEF1D Monoclonal Antibody

单克隆抗体

产品货号: MA01957

第三版

免疫原:Fusion protein of EEF1D

亚型:IgG1

克隆:440

描述:EEF1D,also named as EF1D and EF 1 delta,belongs to the EF-1-beta/EF-1-delta family. It is a subunit of the elongation factor-1 complex,which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. EF-1-beta and EF-1-delta stimulate the exchange of GDP bound to EF-1-alpha to GTP. EEF1D is phosphorylated upon DNA damage,probably by ATM or ATR. This antibody is a mouse monoclonal antibody raised against residues near the C terminus of human EEF1D.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF