



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FAM3C Monoclonal Antibody

单克隆抗体

产品货号: MA01965

第三版

免疫原:Fusion protein of FAM3C

亚型:IgG1

克隆:987

描述:FAM3C,also known as ILEI,is a secreted protein belonging to the cytokine-like protein family. FAM3C is a signaling cytokine with a GG domain,hydrophobic leader sequence and an N-terminal signal peptide. A change in expression of this protein has been noted in pancreatic cancer-derived cells. FAM3C has been involved in the epithelial to mesenchymal transition,and retinal laminar formation processes in vertebrates. Also,FAM3C may play an important role in inner ear development and thus be involved in cell differentiation and proliferation during inner ear embryogenesis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF