



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AHSG Monoclonal Antibody

单克隆抗体

产品货号: MA01967

第三版

免疫原:Fusion protein of Fetuin-A

亚型:IgG2a

克隆:247

描述:Fetuin-A (alpha-2-HS-glycoprotein,AHSG),a liver borne plasma protein,contributes to the prevention of soft tissue calcification,modulates inflammation,reduces insulin sensitivity and fosters weight gain following high fat diet or ageing. It is involved in several functions,such as endocytosis,brain development and the formation of bone tissue. The protein is commonly present in the cortical plate of the immature cerebral cortex and bone marrow hemopoietic matrix,and it has therefore been postulated that it participates in the development of the tissues.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF