

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FHIT Monoclonal Antibody

单克隆抗体

产品货号: MA01968

第三版

免疫原:Fusion protein of FHIT

亚型:IgG1

克隆:308

描述:The candidate tumor suppressor gene,FHIT(Fragile histidine triad protein),encompasses the common human chromosomal fragile site at 3p14.2,the hereditary renal cancer translocation breakpoint,and cancer cell homozygous deletions. Fhit hydrolyzes dinucleotide 5',5'''-P1,P3-triphosphate in vitro and mutation of a central histidine abolishes hydrolase activity. FHIT has a major role in regulating beta-catenin-mediated gene transcription.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF