

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse FN1 Monoclonal Antibody

单克隆抗体

产品货号: MA01972

第三版

免疫原:Fusion protein of Fibronectin

亚型:IgG1

克隆:46

描述:Fibronectin 1 (FN1) is a high molecular weight glycoprotein which exists in both a soluble form in plasma (plasma FN1) and other body fluids and an insoluble form in the extracellular matrix (cellular FN1). Plasma FN1 (dimeric form) is secreted by hepatocytes. Cellular FN (dimeric or cross-linked multimeric forms),made by fibroblasts,epithelial and other cell types,is deposited as fibrils in the extracellular matrix. FN1 binds to cell surfaces through integrins and to various compounds including collagen,fibrin and heparin. It is involved in cell adhesion and migration processes including embryogenesis,wound healing,hemostasis,host defense,and metastasis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF