

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Fibulin 5 Monoclonal Antibody

单克隆抗体

产品货号: MA01973

第三版

免疫原:Fusion protein of Fibulin 5

亚型:IgG1

克隆:937

描述:Fibulin-5,also known as EVEC and DANCE,is a extracellular matrix glycoprotein,secreted by various cell types such as vascular smooth muscle cells, fibroblasts and endothelial cells. Structurally,Fibulin-5 contains an Arg-Gly-Asp (RGD) motif and calcium-binding EGF-like domains. It promotes adhesion of endothelial cells through interaction of integrins and the RGD motif. FBLN5 mutations have been described in patients with age-related macular degeneration,as well as being involved in Charcot-Marie-Tooth neuropathies.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF