

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FOLR1 Monoclonal Antibody

单克隆抗体

产品货号: MA01977

第三版

免疫原:Fusion protein of FOLR1

亚型:IgG1

克隆:783

**描述:**Folate receptor 1 (FOLR1),also known as folate receptor alpha or adult folate-binding protein (FBP),is a 38-kDa glycoprotein belonging to the folate receptor family . The receptor binds to folate and reduced folic acid derivatives and mediates delivery of 5-methyltetrahydrofolate to the interior of cells. FOLR1 is a secreted protein that either anchors to membranes via a glycosyl-phosphatidylinositol linkage or exists in a soluble form. FOLR1 expression is often limited to the apical surfaces of epithelium in the lung,kidney and choroid plexus but is differentially overexpressed in a variety of solid tumors such as ovarian cancer,non-small cell lung cancer,breast cancer,kidney cancer and high-grade osteosarcoma .

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF