



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Gastrin Monoclonal Antibody

单克隆抗体

产品货号: MA01981

第三版

免疫原:Fusion protein of Gastrin

亚型:IgG2a

克隆:404

描述:Gastrin is a hormone produced primarily by G-cells in the stomach,where it functions to stimulate acid secretion by gastric parietal cells. Gastrin is a promiscuous hormone. Thus,the gastrin gene is expressed in common cancers such as bronchogenic,colorectal,ovarian and pancreatic carcinomas.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF