

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat GFAP Monoclonal Antibody

单克隆抗体

产品货号: MA01985

---

第三版

免疫原:Fusion protein of GFAP

亚型:IgG2a

克隆:650

描述:GFAP (Glial fibrillary acidic protein ),an intermediate-filament (IF) protein ,is specifically expressed in cells of astroglial lineage and is widely used to mark the astroglia in the brain. It is also used as a marker for intracranial and intraspinal tumors arising from astrocytes.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF