

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GM2A Monoclonal Antibody

单克隆抗体

产品货号: MA01987

第三版

免疫原:Fusion protein of GM2A

亚型:IgG2a

克隆:347

描述:GM2A,also named as GM2-AP and SAP-3,is the large binding pocket which can accommodate several single chain phospholipids and fatty acids. GM2A also exhibits some calcium-independent phospholipase activity. It plays a key role in the degradation of ganglioside GM2 (GM2). GM2A stimulates only the breakdown of ganglioside GM2 and glycolipid GA2 by beta-hexosaminidase A.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF