



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GOSR2 Monoclonal Antibody

单克隆抗体

产品货号: MA01988

第三版

免疫原:Fusion protein of GOSR2/Membrin

亚型:IgG2a

克隆:99

描述:Golgi SNAP receptor complex member 2 (GOSR2),also known as GS27 or Membrin,is a member of the soluble NSF attachment protein receptor (SNARE) family of vesicle docking proteins. GOSR2 is localized in Golgi apparatus. It is involved in transport of proteins from the cis/medial-Golgi to the trans-Golgi network.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF