



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ADH5 Monoclonal Antibody

单克隆抗体

产品货号: MA01990

第三版

免疫原:Fusion protein of GSNOR/ADH5

亚型:IgG1

克隆:34

描述:ADH5,also named as ADHX,FDH FALDH GSH-FDH ADH3 and GSNOR,belongs to the zinc-containing alcohol dehydrogenase family and Class-III subfamily. It is remarkably ineffective in oxidizing ethanol,but it readily catalyzes the oxidation of long-chain primary alcohols and the oxidation of S-(hydroxymethyl) glutathione. ADH5 mediates multiple cardiovascular functions. It plays in regulating heterocellular communication in the artery wall. . ADH5 immunostaining is distributed in both the nucleus and cytoplasm of the retinal ganglion cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF