



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat HDGF Monoclonal Antibody

单克隆抗体

产品货号: MA01993

第三版

免疫原:Fusion protein of HDGF

亚型:IgG1

克隆:1230

描述:HDGF,also named as HMG1L2,belongs to the HDGF family. It is a heparin-binding protein,with mitogenic activity for fibroblasts. HDGF acts as a transcriptional repressor. It plays important roles in the development and progression of cancers. HDGF may have a synergistic effect with VEGF during cancer development and Progression.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF