

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat HSP90 alpha Monoclonal Antibody

单克隆抗体

产品货号: MA01995

第三版

免疫原:Fusion protein of HSP90

亚型:IgG2a

克隆:1230

**描述:**HSP90,encoded by HSP90AA1,is a constitutively and ubiquitously expressed molecular chaperone that is crucial for the stability and function of many proteins. HSP90 provides chaperoning activity for client proteins; many of them are members of oncogenic pathways,indicating its implication in tumor malignancy. HSP90 mainly resides in the cytosol,while it can also be released to the extracellular space. Secreted Hsp90 is a C-terminal truncated form. It has been reported that the level of plasma Hsp90 is positively correlated with tumor malignancy in clinical cancer patients,and can be a promising diagnostic marker for tumor malignancy in clinical application.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF