

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat HSPA2 Monoclonal Antibody

单克隆抗体

产品货号: MA01996

第三版

免疫原:Fusion protein of HSPA2

亚型:IgG1

克隆:101

描述:Human HSPA2 is a member of the HSPA (HSP70) family of heat-shock proteins,encoded by the gene originally described as testis-specific. Recently,it has been reported that HSPA2 can be also expressed in human somatic tissues in a cell-type specific manner. HSPA2 is a cytosol/nuclear protein able to translocate between cytoplasm and nucleus. HSPA2 is involved in sperm maturity,function and fertility. Aberrant expression of HSPA2 in testes has been connected with male infertility. Recently,HSPA2 has attracted increased interest due to its possible involvement in carcinogenesis of non-testicular tissues. This antibody well recognized the endogenous HSPA4 protein in testis and multiple cell lines. (21373891)

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF