



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PTPRN Monoclonal Antibody

单克隆抗体

产品货号: MA01997

第三版

免疫原:Fusion protein of IA-2/PTPRN

亚型:IgG3

克隆:94

描述:IA-2/PTPRN is a member of the protein tyrosine phosphatase (PTP) family, which contains a transmembrane region, an intracellular PTP-like domain and an extracellular N-terminus. Experiments found that IA-2/PTPRN localizes to secretory granules and is exclusively expressed in neuroendocrine cells (including pancreatic islets cell). IA-2/PTPRN was found to be a major autoantigen in insulin-dependent diabetes mellitus and the detection of autoantibodies against IA-2/PTPRN is commonly used as a diabetes diagnosis marker.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF