

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IDH1 Monoclonal Antibody

单克隆抗体

产品货号: MA01998

第三版

免疫原:Fusion protein of IDH1

亚型:IgG1

克隆:606

描述:IDH1,also named as PICD and IDP,belongs to the isocitrate and isopropylmalate dehydrogenases family. It is a common feature of a major subset of primary human brain cancers. It can form a homodimer.IDH1 mutation is always heterozygotic and IDH1 functions as a dimer,theoretically there will be 25% each wild type and mutant homo-dimers and 50% hetero-dimers present in the tumor cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF