

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IDS Monoclonal Antibody

单克隆抗体

产品货号: MA01999

第三版

免疫原:Fusion protein of Iduronate 2 sulfatase

亚型:IgG1

克隆:3210

**描述:**Iduronate 2-sulfatase (IDS) is required for the lysosomal degradation of heparan sulfate and dermatan sulfate. Mutations causing IDS deficiency in humans result in the lysosomal storage of these glycosaminoglycans and Hunter syndrome,an X chromosome-linked disease. An internal proteolytic cleavage occurs to produce the mature IDS present in human liver shown to contain a 42 kDa polypeptide N-terminal to a 14 kDa polypeptide..

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF