

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IFITM2 Monoclonal Antibody

单克隆抗体

产品货号: MA02001

第三版

免疫原:Fusion protein

亚型:IgG1

克隆:12

描述:IFITM2,also named as 1-8D,belongs to the CD225 family. It is an IFN-induced antiviral protein that mediates cellular innate immunity to at least three major human pathogens,namely influenza A H1N1 virus,West Nile virus (WNV),and dengue virus (WNV),by inhibiting the early step(s) of replication. IFITM2 induces cell cycle arrest and mediates apoptosis by caspase activation and in p53-independent manner. It is overexpressed in colon carcinoma. IFITM2 is a novel pro-apoptotic gene that will provide new insights into the regulated cellular pathways to death. This antibody is specific to IFTIM2.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF