



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse IGBP1 Monoclonal Antibody

单克隆抗体

产品货号: MA02003

第三版

免疫原:Fusion protein of IGBP1

亚型:IgG1

克隆:4789

描述:Immunoglobulin-binding protein 1 (IGBP1),also known as ALPHA-4,was originally identified as a signal transduction molecule coprecipitating with MB1 (CD79a) protein of the B-cell antigen receptor (BCR) complex. IGBP1 has been characterized as an essential regulator of the catalytic activity of protein phosphatase 2A (PP2A),which is a ubiquitous phosphatase with broad substrate specificity and diverse cellular functions.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF