



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse IL1R2 Monoclonal Antibody

单克隆抗体

产品货号: MA02007

第三版

免疫原:Fusion protein of IL1R2

亚型:IgG2a

克隆:123

**描述:**Interleukin-1 receptor type 2 (IL1R2,also known as CD121b) is a cytokine receptor that belongs to the interleukin 1 receptor family. It exists in both a membrane bound and soluble form. IL1R2 binds interleukin alpha (IL1A),interleukin beta (IL1B),and interleukin 1 receptor,type 1(IL1R1/IL1RA),and acts as a decoy receptor that inhibits the activity of its ligands. IL1R2 serves as a negative regulator of IL-1 signaling,thereby suppressing inflammatory response.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF