

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL20RB Monoclonal Antibody

单克隆抗体

产品货号: MA02008

第三版

免疫原:Fusion protein of IL20RB

亚型:IgG2b

克隆:954

**描述:**Interleukin-20 receptor subunit beta (IL20RB) is a single-pass type I membrane protein of the type II cytokine receptor family. The interleukin-20-receptor I complex (IL-20-RI) is composed of two chains,IL20RA and IL20RB. IL-20-RI is a receptor for IL-19,IL-20 and IL-24. IL20RB can also form a heterodimer with IL22RA1. The IL22RA1/IL20RB dimer is a receptor for IL20 and IL24.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF