



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL22RA2 Monoclonal Antibody

单克隆抗体

产品货号: MA02009

第三版

免疫原:Fusion protein of IL22RA2

亚型:IgG2a

克隆:456

**描述:**Interleukin-22 receptor subunit alpha-2 (IL22RA2) is a soluble member of the type II cytokine receptor family. Through binding to IL22,IL22RA2 prevents the interaction of IL22 with the functional IL22-R complex and blocks the activity of IL22 (in vitro). IL22RA2 may play an important role as an IL22 antagonist in the regulation of inflammatory responses. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF