

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL28A Monoclonal Antibody

单克隆抗体

产品货号: MA02010

第三版

免疫原:Fusion protein of IL-28A

亚型:IgG1

克隆:9870

描述:IL28A,also named as IFNL2 and ZCYT020,is a cytokine with immunomodulatory activity. It up-regulates MHC class I antigen expression. IL28A is a ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IFNLR1. The ligand/receptor complex seems to signal through the Jak-STAT pathway. It plays a significant role in the antiviral immune defense in the intestinal epithelium.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF