

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL7 Monoclonal Antibody

单克隆抗体

产品货号: MA02013

第三版

免疫原:Fusion protein of IL-7

亚型:IgG2b

克隆:227

描述:Interleukin-7 (IL7) is a cytokine important for B and T cell development. IL7 plays an active role in the development,survival,maintaining and restoring homeostasis of mature T lymphocytes. It is also a key regulator of the commitment,survival,proliferation and maturation of B cells during development. Furthermore,IL7 can improve the antiviral function and expansion of natural killer (NK) cells and regulate the development and differentiation of dendritic cells. IL7 has also been reported as a regulator of the development of central nervous system and myogenesis and skeletal muscle cell migration.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF