

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat KCNN4 Monoclonal Antibody

单克隆抗体

产品货号: MA02018

---

第三版

免疫原:Fusion protein of KCNN4

亚型:IgG2a

克隆:

描述:The intermediate conductance calcium-activated potassium channel KCNN4,also known as KCa3.1,IK1,or SK4,is a member of the potassium channel KCNN family. KCNN4 is expressed in T cells,mast cells,macrophages,erythrocytes,vascular smooth muscle cells,airway smooth muscle cells,intestinal and airway epithelial cells . Functionally,KCNN4 plays a key role in calcium-dependent cell functions,such as proliferation,activation,and migration,in a broad range of cell types .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF