

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat KIFAP3 Monoclonal Antibody

单克隆抗体

产品货号: MA02019

第三版

免疫原:Fusion protein of KIFAP3

亚型:IgG2a

克隆:456

描述:KIFAP3,also known as KIF3AP or KAP3,is a novel KIF3A/3B-associated protein. It binds to the tail domain of KIF3A/3B and may play a role in regulating the binding of KIF3A/3B to cargoes. Recently it has been reported that mutation within the KIFAP3 gene is associated with decreased KIFAP3 expression and increased survival in sporadic ALS,which makes it a potential target for ALS therapy. KIFAP3 has also been identified as the target of mir-130a.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF