



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human KNG1 Monoclonal Antibody

单克隆抗体

产品货号: MA02020

第三版

免疫原:Fusion protein of Kininogen 1

亚型:IgG1

克隆:834

**描述:**Kininogens are inhibitors of thiol proteases. Kininogen 1 plays important role in Kinin-kallikrein system. This gene is translated into High-molecular weight kininogen (HMWK) and low-molecular weight kininogen (LMWK) after alternative splicing. HMWK is produced by the liver together with prekallikrein. It acts mainly as a cofactor on coagulation and inflammation,and has no intrinsic catalytic activity. LMWK is produced locally by numerous tissues,and secreted together with tissue kallikrein.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF