

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human KLF15 Monoclonal Antibody

单克隆抗体

产品货号: MA02021

第三版

免疫原:Fusion protein of KLF15

亚型:IgM

克隆:9990

描述:KLF15,a kidney-enriched Kr ü ppe-like zinc finger-containing transcription factors,is a transcriptional activator that binds the CLCNKA promoter. The nuclear protein KLF15,is expressed primarily in liver,heart,skeletal muscle and kidney tissues but is not detected in lymphoid tissues or bone marrow. It is an important regulator of GLUT4 in both adipose and muscle tissues..

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF