

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat KMO Monoclonal Antibody

单克隆抗体

产品货号: MA02023

第三版

免疫原:Fusion protein of Kynurenine 3 monooxygenase

亚型:IgG1

克隆:6548

描述:KMO(Kynurenine 3-monooxygenase) is an NADPH-dependent flavin monooxygenase,catalysing the hydroxylation of the L-tryptophan metabolite L-kynurenine to form L-3-hydroxykynurenine. KMO is a membrane protein located on the outer membrane of mitochondria. Tissue distribution studies have revealed that,in rats,highest enzyme activity is found in kidney and liver,with brain having the least activity in comparison to peripheral organs.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF