

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LCK Monoclonal Antibody

单克隆抗体

产品货号: MA02024

第三版

免疫原:Fusion protein of LCK

亚型:IgG1

克隆:

描述:LCK(Lymphocyte cell-specific protein-tyrosine kinase) is also named as LSK,p56-LCK and a nonreceptor protein-tyrosine kinase of the SRC oncogene family that is involved in transduction of T-cell receptor-mediated activation. It is implicated in the control of lymphocyte growth by virtue of its overexpression in some lymphoid malignancies and its transforming activity in heterologous systems. It has 3 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF