

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LGALS3BP Monoclonal Antibody

单克隆抗体

产品货号: MA02025

第三版

免疫原:Fusion protein of LGALS3BP

亚型:IgG2b

克隆:1457

**描述:**The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. LGALS3BP is a widely expressed,secreted galectin. It has been found elevated in the serum of patients with various types of breast,lung,colorectal,ovary,and endometrial cancer,and in those infected by the human immunodeficiency virus (HIV). High expression levels of LGALS3BP are associated with a shorter survival,the occurrence of metastasis or a reduced response to chemotherapy. The native protein binds specifically to a human macrophage-associated lectin known as Mac-2 and also binds galectin 1. It appears to be implicated in immune response associated with natural killer (NK) and lymphokine-activated killer (LAK) cell cytotoxicity. Immunoblotting showed that LGALS3BP is likely contained in breast milk,serum,semen,saliva,urine,and tears.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF