



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LIN28A Monoclonal Antibody

单克隆抗体

产品货号: MA02026

第三版

免疫原:Fusion protein of LIN28

亚型:IgG1

克隆:8700

描述:LIN28 is one of the four key human factors (OCT4,SOX2,NANOG and LIN28) used to reprogram human fibroblasts to an embryonic stem (ES) cell-like state known as the induced pluripotent stem (Ips) cell. Lin28 is a marker of undifferentiated human embryonic stem cells and a cytoplasmic Mrna-binding protein that binds to and enhances the translation of the IGF2 Mrna. LIN28 has also been shown to bind to the let-7 pre-miRNA and block production of the mature let-7 microRNA in mouse embryonic stem cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF