



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human LIPH Monoclonal Antibody

单克隆抗体

产品货号: MA02027

第三版

免疫原:Fusion protein of LIPH

亚型:IgG2a

克隆:1246

描述:LIPH(Lipase member H) is also named as LPDLR,MPAPLA1,PLA1B. It is a phosphatidic acid-selective phospholipase A1 (PLA1) that produces 2-acyl lysophosphatidic acid (LPA) and belongs to the AB hydrolase superfamily and Lipase family.LIPH has a 12-residue lid region,which likely covers a catalytic pocket,and 4 potential N-linked glycosylation sites. Defects in LIPH are the cause of hypotrichosis type 7 (HYPT7).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF