

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat LONP1 Monoclonal Antibody

单克隆抗体

产品货号: MA02028

第三版

免疫原:Fusion protein of LONP1

亚型:IgG1

克隆:3200

描述:LONP1(Lon protease homolog,mitochondrial) is also named as LONP,LONHS,HLON,LON,PRSS15,PIM1,MGC1498 and belongs to the peptidase S16 family. It seems to play a major role in the elimination of oxidatively modified proteins in the mitochondrial matrix. LONP1,also a nuclearly encoded and mitochondrially located stress-responsive protease,is involved in heme-mediated ALAS-1 turnover. It recognizes specific surface determinants or folds,initiates proteolysis at solvent-accessible sites,and generates unfolded polypeptides that are then processively degraded.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF