



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse LYN Monoclonal Antibody

单克隆抗体

产品货号: MA02030

第三版

免疫原:Fusion protein of LYN

亚型:IgG2a

克隆:235

描述:LYN,also named as JTK8,p53Lyn and p56Lyn,belongs to the protein kinase superfamily,Tyr protein kinase family and SRC subfamily. It is a non-receptor tyrosine-protein kinase that transmits signals from cell surface receptors and plays an important role in the regulation of innate and adaptive immune responses,hematopoiesis,responses to growth factors and cytokines,integrin signaling,but also responses to DNA damage and genotoxic agents. LYN is a component of the mTOR (mammalian target of rapamycin) pathway.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF