

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MAD2L1 Monoclonal Antibody

单克隆抗体

产品货号: MA02031

第三版

免疫原:Fusion protein of MAD2L1

亚型:IgG2b

克隆:432

描述:MAD2 mitotic arrest deficient-like 1 (yeast) (MAD2L1,synonyms: MAD2,HSMAD2) is a component of the mitotic spindle assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly aligned at the metaphase plate. MAD2L1 is localized to human chromosome band 5q23.3 by in situ hybridization.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF