



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SCGB2A2 Monoclonal Antibody

单克隆抗体

产品货号: MA02032

第三版

免疫原:Fusion protein of Mammaglobin

亚型:IgG1

克隆:

描述:Mammaglobin A,also known as SCGB2A2,is a member of the secretoglobin superfamily. Mammaglobin A is a breast cancer-associated antigen almost exclusively over-expressed in primary and metastatic human breast cancers,making it a specific molecular marker and a potential therapeutic target for breast cancer. This monoclonal antibody is raised against full-length of 93-amino acid human mammaglobin A.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF