

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat MCRS1 Monoclonal Antibody

单克隆抗体

产品货号: MA02034

---

第三版

免疫原:Fusion protein of MCRS1

亚型:IgG2a

克隆:332

**描述:**MCRS1,also named INO80Q,MSP58 or p78,is a nucleolar protein expressed in the very early S phase. It modulates the transcription repressor activity of DAXX by recruiting it to the nucleolus. MCRS1 is also present on polyribosomes of synaptoneurosome as a novel RNA-binding protein to activate the rRNA transcription. The forkhead-associated domain of MCRS1 is involved in association with centrosomal protein Nde1,comprising a component of the centrosome and taking an essential role in viability. Its isoform MCRS2 might be a linker between telomere maintenance and cell-cycle regulation by interacting with LPTS/PinX1,a telomerase-inhibitory protein. MCRS1 is detected by Western blot as a highly modified protein at 78 kDa and minor bands at about 62 and 55 kDa.

**配方:**

**浓度:**None

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF