

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MLPH Monoclonal Antibody

单克隆抗体

产品货号: MA02035

第三版

免疫原:Fusion protein of Melanophilin

亚型:IgG1

克隆:72

描述:Melanophilin (MLPH) is a member of the exophilin subfamily of Rab effector proteins. MLPH is a modular protein,forming a ternary complex with the small Ras-related GTPase Rab27A in its GTP-bound form and the motor protein myosin Va. It regulates retention of melanosomes at the peripheral actin cytoskeleton of melanocytes,a process essential for normal mammalian pigmentation. Defects in MLPH are a cause of Griscelli syndrome type 3 (GS3),which is characterized by a silver-gray hair color and abnormal pigment distribution in the hair shaft.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF