



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD10 Monoclonal Antibody

单克隆抗体

产品货号: MA02038

第三版

免疫原:Fusion protein of MME/CD10

亚型:IgG2b

克隆:9870

描述:MME is also names as neprilysin,CALLA,NEP,SFE,CD10 and belongs to the peptidase M13 family.It is a glycoprotein that is particularly abundant in kidney,where it is present on the brush border of proximal tubules and on glomerular epithelium.CD10 has been considered a useful marker in the diagnosis of renal carcinomas,because of its expression in clear cell and papillary renal cell carcinomas and its absence in chromophobe renal cell carcinomas.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF